



ZENEX
MULTI
I-SYSTEM

IZEN Implant,

bright smile

Certification



History & Global Network



2016. 03

Foundation of IZEN Implant

2016. 06

Acquisition of manufacturing licenses for medical devices
Report medical device items (Screw driver for dental implant surgery)

2016. 12

Patent registration (Driver holder and driver apparatus for dental engineering, No. 10-1684959)

2017. 03

Patent registration (Driver holder for dental apparatus, No. 10-1714533, No. 10-1714541)

2018. 08

Venture Business Certificate (No. 20180110278)

2018. 09

Foundation of Corporate Affiliated Research Institute (No. 2018114470)

2019. 06

Acquisition of small businesses technology development project (Development of Dental implant abutments and scan body for dental digital system)

2019.12

Patent registration (Abutment and dental implant, No. 10-2062575)



IZEN

I M P L A N T

Your companion armed with active
and flexible ideas for your satisfaction.

2017.09

Registration of medical devices
(Apparatus for dental implant
surgery)

2018.03

Certification of medical
device items compliance with
manufacturing and quality
control standards (KGMP)

2018.04

Medical Device Item Certification
(Dental implant abutments)

2018.07

Acquisition of ISO 9001 : 2015

2020.01

Medical Device Item
Certification
(Orthodontic Screw)

2020.12

Medical Device Item
Certification
(Dental implant abutments)

2021.01

Patent registration
(Implant Lab Analog,
No. 10-2209274)

2021.05

Acquisition of ISO13485:
2016 & CE MDD

MANUFACTURE

Manufacturing



Materials

Fixture and Abutment of ZENEX MULTI & ZENEX PLUS Implant System are produced by titanium which is certified by ASTM.



Manufacture by CNC Machine

In order to produce Fixture and Abutment of ZENEX MULTI & ZENEX PLUS Implant System, the raw materials are manufactured by CNC Machine.




Manufacture of Half-Finished product

Fixture and Abutment of ZENEX MULTI & ZENEX PLUS Implant System are manufactured as a Half-Finished product based on a design of each product.



Sand Blasting

In case of Fixture, the roughness of surface is increased physically by blasting Alumina powder for osseointegration through getting a wide SLA surface.



Etching

In case of Fixture, the roughness of surface is widened by chemical corrosion. After a SLA surface treatment, the surface inspection is done whether any residual acid is on the fixture.



TiN Coating

For classifying product and aesthetic impression, the electronic-chemical methods are used artificially for TiN coating as a gold color. Applied for some kinds of Abutments.



Cleaning

Final cleaning for a perfect removal of pollutant is proceeded in the Clean Room. After final cleaning, the inspection will be taken place.



Packing

After cleaning products, the packing is done in the Clean Room. The packing is proceeded by automatic system without any foreign pollution.



Sterilization

For some products, it is sterilized by the short-wavelength of gamma radiation. It is possible to check the sterilization through a color of sticker on the package changed by a gamma radiation. Applied product: Fixture, Cover Screw, Healing Abutment & Multi Healing Cap.

Quality Management



Dimension Management

After manufacturing product (fixture), Dimension of product is tested by non-contact visual microscope to measure precise dimension.



Stability Test

Our products are tested by fixation force and screw loosening for providing a stable products.




Fatigue Test

For providing the strength of product and a long-term safety, a fatigue test is done by measuring a mechanical force. Based on ISO14801:2016, the test is taken place.



Surface Management

In order to maintain a stable roughness of SLA Implant surface, the roughness of surface is tested in each lot.



Connection Management

We observe the fixture and abutment is connected internally without any problem before the shipment.



ZENEX MULTI I-System

COMPOSITION

022 ZENEX MULTI I-System Fixture

025 Cover Screw

026 Healing Abutment

029 PROSTHETIC FLOW CHART I

030 Temporary Abutment

031 Cemented Abutment

036 Angled Abutment

040 FreeMilling Abutment

042 CCM Cast Abutment

043 Pick-Up Impression Coping

045 Transfer Impression Coping

047 Fixture Lab Analog

049 PROSTHETIC FLOW CHART II

050 Multi Straight Abutment

051 Multi Angled Abutment

053 Components of Multi Abutment

059 PROSTHETIC FLOW CHART III

060 Ball Abutment

061 Components of Ball Abutment

064 Kerator Abutment

065 Kerator Angled Abutment

066 Components of Kerator Abutment

070 PROSTHETIC FLOW CHART IV

071 ScanBody


072 Ti Link Abutment

073 Ti Blank Abutment

ZENEX MULTI I-System

- Design features

- Designed for various types of bone



ZENEX MULTI Implant

- Surface features

SEM MAG:30X

SLA surface
(Al₂O₃ powder Sand-blast and Acid Etched)

- SLA surface for presenting ideal surface shape
- Surface Roughness Ra 2.0~3.0 μm
- Safe surface with no residual acid
- Safer than other implants
(Proved by ICP/IC Analysis)

SEM MAG:5000X

ZENEX MULTI Fixture I-System

Internal 11° tapered submerged type implant
 SLA surface implant
 Tapered design to secure initial stability
 Double thread design for securing width of internal connection and easy insertion
 Possible to secure primary stability on apex of fixture by realizing wide apical diameter.
 Recommended torque for implant placement: Less than 40Ncm

For Narrow Fixture (Ø 3.0)

- Enable surgery in bone of narrow width
- Easy to compensate angle in anterior area

For Wide Fixture (Ø 6.0, Ø 7.0)


- Wide fixture with single thread for easy implant placement.
- 6mm length fixture enables to use in lower alveolar bone
- Radial shaped on upper part of fixture enables to be used for wide hole such as extraction socket

* Recommend to use fixture bigger than Ø4.5 mm for single case in posterior area



* Incompatible to use other brands' prosthetic products for Short Fixture (less than 7mm length) in case of I-System

Cover Screw not included in the Fixture package.



022 ZENEX MULTI





R
Ø4.5
(mm)





W
Ø5.0
(mm)





W
Ø5.5
(mm)



W
Ø6.0
(mm)








W
Ø7.0
(mm)






023 ZENEX
MULTI



Special Options

<p>Narrow Ø3.0 (mm)</p> 		15	<p>Mini Ø3.5 (mm)</p> 		15
		ZMSN3015			ZMSM3515

<p>R Ø4.0 (mm)</p> 		15	18	20
		ZMSR4015	ZMSR4018	ZMSR4020

<p>R Ø4.5 (mm)</p> 		15	18	20
		ZMSR4515	ZMSR4518	ZMSR4520

<p>W Ø5.0 (mm)</p> 		6	15
		ZMSW5006	ZMSW5015










<p>W Ø5.5 (mm)</p> 		6	15
		ZMSW5506	ZMSW5515

024 ZENEX
MULTI

Cover Screw

Select appropriate fixture height upon depth of implant placement.
Select specification fits for fixture connection.



Tighten with 1.2 Hex Driver by hand
Recommended tightening torque: 5~8Ncm

	H (mm)	0.4	1.4	2.0
Narrow				
		CSSN30S	CSSN30M	CSSN30L
Mini				
		CSSM30S	CSSM30M	CSSM30L
R/W				
		CSSR40S	CSSR40M	CSSR40L

Healing Abutment

Use Healing Abutment fits for the diameter of abutment.
 Use specification fits for fixture connection.





Tighten with 1.2 Hex Driver by hand
 Recommended tightening torque: 5-8Ncm






Narrow

H (mm)	3.0	4.0	5.0	7.0	9.0
Ø4.3	HATM403	HATM404	HATM405	HATM407	HATM409
Ø4.8	HATM453	HATM454	HATM455	HATM457	HATM459

M/R/W

H (mm)	2.0	3.0	4.0	5.0
D (mm)				
Ø4.3	HASR402	HASR403	HASR404	HASR405
Ø4.8	HASR452	HASR453	HASR454	HASR455
Ø5.5	HASR522	HASR523	HASR524	HASR525
Ø6.0	HASR572	HASR573	HASR574	HASR575
Ø6.8	HASR652	HASR653	HASR654	HASR655
Ø8.0		HASR803	HASR804	HASR805
Ø9.0		HASR903	HASR904	HASR905

M/R/W


H (mm)	6.0	7.0	9.0
D (mm)			
Ø4.3	HASR406	HASR407	HASR409
Ø4.8	HASR456	HASR457	HASR459
Ø5.5	HASR526	HASR527	HASR529
Ø6.0	HASR576	HASR577	HASR579
Ø6.8	HASR656	HASR657	HASR659

028 ZENEX
MULTI

IZEN

I M P L A N T

Prosthetic Flow Chart I




Temporary Abutment

Abutment for manufacturing Screw-retained type temporary prosthesis
 Select specification fits for fixture connection.
 Fixture Level Impression


Tighten with 1.2 Hex Driver

Recommended tightening torque: 20Ncm


To order as a set as 'Abutment + Abutment Screw'
 : Product Code + S (ex: TASRH4510S)



Narrow
 DØ4.0



M/R/W
 DØ4.5



030 ZENEX MULTI

Cemented Abutment

Abutment for manufacturing Cement/Combination-retained type prosthesis

Select specification fits for fixture connection.


Customized by grinding (need to be maintained at least 3.0mm of Abutment Length above Fixture Platform)

Fixture Level Impression

Tighten with 1.2 Hex Driver

Recommended tightening torque: 30Ncm

To order as a set as 'Abutment + Abutment Screw'
: Product Code + S (ex: CASRH45104S)



Narrow DØ4.0

		G/H (mm)	1.0	2.0	3.0	4.0
Hex	P/H (mm)					
	5.5	CATMH40105	CATMH40205	CATMH40305	CATMH40405	
	7.0	CATMH40107	CATMH40207	CATMH40307	CATMH40407	
Non-Hex	5.5	CATMN40105	CATMN40205	CATMN40305	CATMN40405	
	7.0	CATMN40107	CATMN40207	CATMN40307	CATMN40407	

Narrow DØ4.0

		G/H (mm)	5.0	6.0	7.0
Hex	P/H (mm)				
	5.5	CATMH40505	CATMH40605	CATMH40705	
	7.0	CATMH40507	CATMH40607	CATMH40707	
Non-Hex	5.5	CATMN40505	CATMN40605	CATMN40705	
	7.0	CATMN40507	CATMN40607	CATMN40707	





Narrow
DØ4.5

		G/H (mm)	1.0	2.0	3.0	4.0
		P/H (mm)				
Hex	5.5	CATMH45105	CATMH45205	CATMH45305	CATMH45405	
	7.0	CATMH45107	CATMH45207	CATMH45307	CATMH45407	
Non-Hex	5.5	CATMN45105	CATMN45205	CATMN45305	CATMN45405	
	7.0	CATMN45107	CATMN45207	CATMN45307	CATMN45407	

Narrow
DØ4.5

		G/H (mm)	5.0	6.0	7.0
		P/H (mm)			
Hex	5.5	CATMH45505	CATMH45605	CATMH45705	
	7.0	CATMH45507	CATMH45607	CATMH45707	
Non-Hex	5.5	CATMN45505	CATMN45605	CATMN45705	
	7.0	CATMN45507	CATMN45607	CATMN45707	





M/R/W
DØ4.5

		G/H (mm)	1.0	2.0	3.0	4.0
		P/H (mm)				
Hex	4.0	CASRH45104	CASRH45204	CASRH45304	CASRH45404	
	5.5	CASRH45105	CASRH45205	CASRH45305	CASRH45405	
	7.0	CASRH45107	CASRH45207	CASRH45307	CASRH45407	
Non-Hex	4.0	CASRN45104	CASRN45204	CASRN45304	CASRN45404	
	5.5	CASRN45105	CASRN45205	CASRN45305	CASRN45405	
	7.0	CASRN45107	CASRN45207	CASRN45307	CASRN45407	

M/R/W
DØ4.5

		G/H (mm)	5.0	6.0	7.0
		P/H (mm)			
Hex	4.0	CASRH45504	CASRH45604	CASRH45704	
	5.5	CASRH45505	CASRH45605	CASRH45705	
	7.0	CASRH45507	CASRH45607	CASRH45707	
Non-Hex	4.0	CASRN45504	CASRN45604	CASRN45704	
	5.5	CASRN45505	CASRN45605	CASRN45705	
	7.0	CASRN45507	CASRN45607	CASRN45707	





M/R/W
DØ5.2

		G/H (mm)	1.0	2.0	3.0	4.0
		P/H (mm)				
Hex	4.0	CASRH52104	CASRH52204	CASRH52304	CASRH52404	
	5.5	CASRH52105	CASRH52205	CASRH52305	CASRH52405	
	7.0	CASRH52107	CASRH52207	CASRH52307	CASRH52407	
Non-Hex	4.0	CASRN52104	CASRN52204	CASRN52304	CASRN52404	
	5.5	CASRN52105	CASRN52205	CASRN52305	CASRN52405	
	7.0	CASRN52107	CASRN52207	CASRN52307	CASRN52407	

M/R/W
DØ5.2

		G/H (mm)	5.0	6.0	7.0
		P/H (mm)			
Hex	4.0	CASRH52504	CASRH52604	CASRH52704	
	5.5	CASRH52505	CASRH52605	CASRH52705	
	7.0	CASRH52507	CASRH52607	CASRH52707	
Non-Hex	4.0	CASRN52504	CASRN52604	CASRN52704	
	5.5	CASRN52505	CASRN52605	CASRN52705	
	7.0	CASRN52507	CASRN52607	CASRN52707	

M/R/W
DØ5.7

		G/H (mm)	1.0	2.0	3.0	4.0
		P/H (mm)				
Hex	4.0	CASRH57104	CASRH57204	CASRH57304	CASRH57404	
	5.5	CASRH57105	CASRH57205	CASRH57305	CASRH57405	
	7.0	CASRH57107	CASRH57207	CASRH57307	CASRH57407	
Non-Hex	4.0	CASRN57104	CASRN57204	CASRN57304	CASRN57404	
	5.5	CASRN57105	CASRN57205	CASRN57305	CASRN57405	
	7.0	CASRN57107	CASRN57207	CASRN57307	CASRN57407	

M/R/W
DØ5.7

		G/H (mm)	5.0	6.0	7.0
		P/H (mm)			
Hex	4.0	CASRH57504	CASRH57604	CASRH57704	
	5.5	CASRH57505	CASRH57605	CASRH57705	
	7.0	CASRH57507	CASRH57607	CASRH57707	
Non-Hex	4.0	CASRN57504	CASRN57604	CASRN57704	
	5.5	CASRN57505	CASRN57605	CASRN57705	
	7.0	CASRN57507	CASRN57607	CASRN57707	

M/R/W
DØ6.5

		G/H (mm)	1.0	2.0	3.0	4.0
		P/H (mm)				
Hex	4.0	CASRH65104	CASRH65204	CASRH65304	CASRH65404	
	5.5	CASRH65105	CASRH65205	CASRH65305	CASRH65405	
Non-Hex	4.0	CASRN65104	CASRN65204	CASRN65304	CASRN65404	
	5.5	CASRN65105	CASRN65205	CASRN65305	CASRN65405	

M/R/W
DØ6.5

	G/H (mm)	5.0	6.0	7.0
	P/H (mm)			
Hex	4.0	CASRH65504	CASRH65604	CASRH65704
	5.5	CASRH65505	CASRH65605	CASRH65705
Non-Hex	4.0	CASRN65504	CASRN65604	CASRN65704
	5.5	CASRN65505	CASRN65605	CASRN65705

035
ZENEX
MULTI

Angled Abutment

Abutment for manufacturing Cement/Combination-retained type prosthesis

Various types of Angle (17° for Ø 3.0 Narrow Fixture / 15° & 25° for Mini, Regular and Wide Fixture [Ø 3.5 ~ Ø 7.0])

Select specification fits for fixture connection.

Fixture Level Impression


Can be positioned in 12 directions by selecting A or B type

Tighten with 1.2 Hex Driver

Recommended tightening torque: 30Ncm

To order as a set as 'Abutment + Abutment Screw'

: Product Code + S (ex: AASRA45152S)



Narrow
DØ4.0

G/H (mm)	2.0	4.0
Angle	17°	



A-Type

AATMA40172

AATMA40174

B-Type

AATMB40172

AATMB40174

Non-Hex

AATMN40172

AATMN40174

Narrow
DØ4.5

G/H (mm)	2.0	4.0
Angle	17°	



A-Type

AATMA45172

AATMA45174

B-Type


AATMB45172

AATMB45174

Non-Hex

AATMN45172

AATMN45174



M/R/W
DØ4.5

G/H (mm)	2.0	3.0	4.0
Angle		15°	



A-Type

AASRA45152

AASRA45153

AASRA45154

B-Type

AASRB45152

AASRB45153

AASRB45154

Non-Hex

AASRN45152

AASRN45153

AASRN45154

M/R/W
DØ5.2

G/H (mm)	2.0	3.0	4.0
Angle		15°	



A-Type

AASRA52152

AASRA52153

AASRA52154

B-Type

AASRB52152

AASRB52153

AASRB52154

Non-Hex

AASRN52152

AASRN52153

AASRN52154

M/R/W
DØ5.7

G/H (mm)	2.0	3.0	4.0
Angle	15°		



A-Type

AASRA57152

AASRA57153

AASRA57154

B-Type

AASRB57152

AASRB57153

AASRB57154

Non-Hex

AASRN57152

AASRN57153

AASRN57154

M/R/W
DØ4.5

G/H (mm)	2.0	3.0	4.0
Angle	25°		



A-Type

AASRA45252

AASRA45253

AASRA45254

B-Type

AASRB45252

AASRB45253

AASRB45254

Non-Hex

AASRN45252

AASRN45253

AASRN45254

M/R/W
DØ5.2

G/H (mm)	2.0	3.0	4.0
Angle	25°		



A-Type

AASRA52252

AASRA52253

AASRA52254

B-Type

AASRB52252

AASRB52253

AASRB52254

Non-Hex

AASRN52252

AASRN52253

AASRN52254

M/R/W
DØ5.7

G/H
(mm)
Angle

2.0

3.0

4.0

25°



A-Type

AASRA57252

AASRA57253

AASRA57254

B-Type

AASRB57252

AASRB57253

AASRB57254

Non-Hex

AASRN57252

AASRN57253


AASRN57254

FreeMilling Abutment

Abutment for manufacturing Cement/Combination-retained type prosthesis
 Used when creating free marginal space for Abutment
 Select specification fits for fixture Connection
 Customized by grinding (need to be maintained at least 3.0mm of Abutment Length above Fixture Platform)
 Fixture Level Impression

Tighten with 1.2 Hex Driver
 Recommended tightening torque: 30Ncm
 When using a Narrow Fixture (Ø 3.0), use Mini Type FreeMilling Abutment of T-System.

To order as a set as 'Abutment + Abutment Screw'
 : Product Code + S (ex: FMASRH52208S)



040 ZENEX MULTI

Narrow

DØ4.0
 H: 9.0mm



Hex

FMATMH40159



Non-Hex

FMATMN40159



Hex

FMATMH40309



Non-Hex

FMATMN40309

M/R/W

DØ4.5
 H: 9.0mm



Hex

FMASRH45209



Non-Hex

FMASRN45209



Hex

FMASRH45409



Non-Hex

FMASRN45409

M/R/W
DØ5.2
H: 8.0mm



M/R/W
DØ5.7
H: 8.0mm



M/R/W
DØ6.5
H: 8.0mm




041 ZENEX
MULTI

CCM Cast Abutment


Abutment for manufacturing customized abutment in difficult and complicated cases
 Select specification fits for fixture connection
 Fixture Level Impression
 Casting with non-precious alloy for manufacturing customized prosthesis
 Melting point of CCM : 1,400 ~ 1,550°C

Tighten with 1.2 Hex Driver
 Recommended tightening torque: 30Ncm


To order as a set as 'Abutment + Abutment Screw'
 : Product Code + S (ex: CCASRH4510S)



Narrow DØ4.0



M/R/W DØ4.5




Pick-Up Impression Coping


Components for Fixture Level Impression
 Open Tray Type
 Select specification fits for fixture connection
 Enables accurate impression with design of stably fixed in impression material

Tighten with 1.2 Hex Driver by hand

To order as a set as 'Impression Coping Body + Guide Pin'
 : [For 11mm Length] Product Code + S (ex: PISR40SS)
 [For 15mm Length] Product Code + S (ex: PISR40LS)









Narrow



Mini, Regular & Wide

043 ZENEX MULTI

		L (mm)		11.0		Guide Pin S	
		Hex	Non-Hex				
							
Narrow	D (mm)					PIGPTMS	
	Ø4.0	PITMH40S	PITMN40S				
	Ø4.5	PITMH45S	PITMN45S				
M/R/W	Ø4.0	PISR40S	PISRN40S				
	Ø4.8	PISR45S	PISRN45S				
	Ø5.5	PISR52S	PISRN52S				
	Ø6.0	PISR57S	PISRN57S				
	Ø6.8	PISR65S	PISRN65S			PIGPSRS	


		L (mm)		15.0		Guide Pin L	
		Hex	Non-Hex				
		D (mm)					
Narrow	Ø4.0	PITMH40L	PITMN40L			PIGPTML	
	Ø4.5	PITMH45L	PITMN45L			PIGPTML	
M/R/W	Ø4.0	PISRH40L	PISRN40L				
	Ø4.8	PISRH45L	PISRN45L				
	Ø5.5	PISRH52L	PISRN52L				
	Ø6.0	PISRH57L	PISRN57L				
	Ø6.8	PISRH65L	PISRN65L			PIGPSRL	

Transfer Impression Coping


Components for Fixture Level Impression
 Closed Tray Type
 Select specification fits for fixture connection
 Streamlined Shape: Easy to transfer
 Anti-Rotation Grooves accord with Hex of fixture

Tighten with 1.2 Hex Driver by hand

To order as a set as 'Impression Coping Body + Guide Pin'
 : [For 9mm Length] Product Code + S (ex: TISRH40S)
 [For 13mm Length] Product Code + S (ex: TISRH40LS)









Narrow



Mini, Regular & Wide

045 ZENEX MULTI

	L (mm)	9.0		Guide Pin S	
		Hex	Non-Hex	S	
					
	D (mm)				
Narrow	Ø4.0	TITMH40S	TITMN40S	TIGPTMS	
	Ø4.5	TITMH45S	TITMN45S		
M/R/W	Ø4.0	TISRH40S	TISRN40S	TIGPSRS	
	Ø4.8	TISRH45S	TISRN45S		
	Ø5.5	TISRH52S	TISRN52S		
	Ø6.0	TISRH57S	TISRN57S		
	Ø6.8	TISRH65S	TISRN65S		

		L (mm)		13.0		Guide Pin L	
		Hex	Non-Hex				
		D (mm)					
Narrow	Ø4.0	TITMH40L	TITMN40L			TIGPTML	
	Ø4.5	TITMH45L	TITMN45L			TIGPTML	
M/R/W	Ø4.0	TISRH40L	TISRN40L			TIGPSRL	
	Ø4.8	TISRH45L	TISRN45L			TIGPSRL	
	Ø5.5	TISRH52L	TISRN52L			TIGPSRL	
	Ø6.0	TISRH57L	TISRN57L			TIGPSRL	
	Ø6.8	TISRH65L	TISRN65L			TIGPSRL	

Fixture Lab Analog

Lab Analog for Fixture Level Impression

Select among Fixture Diameter (\varnothing 3.0 / \varnothing 3.5 / \varnothing 4.0 and above)



Narrow
D \varnothing 3.2



FLASN300

Mini
D \varnothing 3.75



FLASM350

R/W
D \varnothing 4.25



FLASR400


047
ZENEX
MULTI

048 ZENEX
MULTI

IZEN

I M P L A N T

Prosthetic Flow Chart II




049
ZENEX
MULTI

Multi Straight Abutment

Abutment for manufacturing screw-retained prosthesis in Multiple Case
 Same platform as Multi Angled Abutment
 Move into internal oral part by using exclusive Abutment Carrier (Code: MSACR48)

Tighten with exclusive driver (Code: MSADSR20)
 Recommended tightening torque: 30Ncm

To order as a set as 'Abutment + Carrier'
 : Product Code + S (ex: MSASR4815S)



Narrow

Mini, Regular & Wide

Narrow DØ4.8

G/H (mm) 1.0 2.0 3.0 4.0 5.0



MSATM4810 MSATM4820 MSATM4830 MSATM4840 MSATM4850

M/R/W DØ4.8

G/H (mm) 1.5 2.5 3.5 4.5



MSASR4815 MSASR4825 MSASR4835 MSASR4845

Multi Angled Abutment

Abutment for manufacturing screw-retained prosthesis in Multiple Case
 Abutment of various angles (17°, 30°) for various angled of implant insertion path
 Same platform as Multi Straight Abutment
 Compensation of fixture placement angle up to 108°
 Connect by using exclusive Abutment Carrier (Code: MAACRMC)

Tighten with 1.2 Hex Driver

Recommended tightening torque: 30Ncm

Multi Angled Abutment Screw (MAASTM20 for Mini & MAASSR23 for Regular and Wide) included

To order as a set as 'Abutment + Screw + Carrier'
 : Product Code + S (ex: MAASRH481725S)



051 ZENEX MULTI

Narrow DØ4.8

G/H (mm)	2.5	3.0	4.0
----------	-----	-----	-----

Angle
17°



MAATMH481725 MAATMH481730 MAATMH481740

Narrow DØ4.8

G/H (mm)	3.5	4.0	5.0
----------	-----	-----	-----

Angle
30°



MAATMH483035 MAATMH483040 MAATMH483050

M/R/W
DØ4.8

G/H
(mm)

2.5

3.5

4.5

5.5

Angle
17°



Hex

MAASRH481725

MAASRH481735

MAASRH481745

MAASRH481755

Non-Hex

MAASRN481725

MAASRN481735

MAASRN481745

MAASRN481755

M/R/W
DØ4.8

G/H
(mm)

3.5

4.5

5.5

6.5

Angle
30°



Hex

MAASRH483035

MAASRH483045

MAASRH483055

MAASRH483065

Non-Hex

MAASRN483035

MAASRN483045

MAASRN483055

MAASRN483065

Multi Temporary Cylinder


Multi Cylinder for making Combination-retained prosthesis from Multi Abutment
For manufacturing Screw-retained temporary Abutment

Tighten with 1.2 Hex Driver

Recommended tightening torque: 20Ncm

Multi Cylinder Screw (MTCSR23) included

To order as a set as 'Cylinder + Multi Cylinder Screw'
: Product Code + S (ex: MTCSR48S)



DØ4.8

MTCSR48

053 ZENEX MULTI

Multi CCM Cast Cylinder


Multi Cylinder for making Screw-retained prosthesis from Multi Abutment
Casting with non-precious alloy for manufacturing customized prosthesis
Melting point of CCM : 1,400 ~ 1,550°C

Tighten with 1.2 Hex Driver

Recommended tightening torque: 20Ncm

Multi Cylinder Screw (MTCSR23) included

To order as a set as 'Cylinder + Multi Cylinder Screw'
: Product Code + S (ex: MCCC48S)



DØ4.8

MCCC48

Multi Plastic Cast Cylinder


Multi Cylinder for making Screw-retained prosthesis from Multi Abutment

Tighten with 1.2 Hex Driver

Recommended tightening torque: 20Ncm

Multi Cylinder Screw (MTCSR23) included

To order as a set as 'Cylinder + Multi Cylinder Screw'
: Product Code + S (ex: MPCCSRN48S)



DØ5.0

MPCCSRN48

Multi Ti Link Cylinder

Multi Cylinder for making Combination-retained prosthesis from Multi Abutment

Use when producing Customized Abutment (Titanium & Zirconia) and Crown by CAD/CAM equipment

Use exclusive implant library of ZENEX MULTI Implant System


Abutment Level Impression

Tighten with 1.2 Hex Driver

Recommended tightening torque: 20Ncm

Multi Cylinder Screw (MTCSR23) included

To order as a set as 'Cylinder + Multi Cylinder Screw'
: Product Code + S (ex: MTLCSRN48S)



DØ4.8

MTLCSRN48


Multi Pick-Up Impression Coping

Components for taking Abutment Level impression for Multi Abutment Open Tray Type

Tighten with 1.2 Hex Driver by hand

Multi Pick-up Impression Coping Guide Pin (MPICSRGP) included

To order as a set as 'Impression Coping body + Guide Pin'
: Product Code + S (ex: MPICSR48S)



DØ4.8


MPICSR48

055
ZENEX
MULTI

Multi Transfer Impression Coping

Components for taking Abutment Level impression for Multi Abutment Closed Tray Type

Tighten with 1.2 Hex Driver by hand



DØ4.8

MTICSR48

Multi Lab Analog

Lab Analog for Multi Abutment



MLASR48

056
ZENEX
MULTI

Multi Healing Cap

Protect Cap for Multi Abutment

Tighten with 1.2 Hex Driver by hand



$D\varnothing 4.8$

MHCSR48

Multi ScanBody

ScanBody for Multi Abutment

ScanBody for manufacturing customized Titanium abutment
Use for Oral scan (Model scan available as well)

Tighten with 1.2 Hex Driver by hand

To order as a set as 'ScanBody + ScanBody Screw'
: Product Code + S (ex: MSBSR4809S)



MSBSR4809

057
ZENEX
MULTI

Multi Straight Abutment Driver

Torque driver for Multi Straight Abutment



MSADSR20

058
ZENEX
MULTI


IZEN

I M P L A N T

Prosthetic Flow Chart III

0 lbs 2.4 lbs 1.2 lbs 4 lbs

Angled Light Angled Extra light Angled




059 ZENEX MULTI

Ball Abutment

Abutment for overdenture using O-ring attachment
 Compensation of mounting angle up to 20°

Tighten with exclusive Ball Abutment Driver (Code: BAD24)
 Recommended tightening torque: 30Ncm



060 ZENEX
 MULTI

G/H (mm)	1.0	2.0	3.0	4.0	5.0	6.0
Narrow DØ3.5						
	BATM3510	BATM3520	BATM3530	BATM3540	BATM3550	BATM3560

G/H (mm)	1.0	2.0	3.0	4.0	5.0	6.0
M/R/W DØ3.5						
	BASR3510	BASR3520	BASR3530	BASR3540	BASR3550	BASR3560

Ball Abutment Cap Set

O-ring attachment for Ball Abutment

Replace O-ring on Ball Abutment Cap

Packing Unit: Ball Abutment Cap + Ball Abutment Lab O-ring



BASRCS

061 ZENEX
MULTI

Ball Abutment Retainer Set

Use in case lack of vertical diameter comparing to Ball Abutment Retainer

Packing Unit: Ball Abutment Retainer + Ball Abutment Lab O-ring



BASRRS

Ball Abutment O-ring Set

O-ring set

Packing Unit: Ball Abutment O-ring 5 EA



BASROS

Ball Abutment Lab Analog

Lab Analog for Ball Abutment



BAALA

Ball Abutment Driver

Torque driver for Ball Abutment




BAD24

063
ZENEX
MULTI

Kerator Abutment

Compensation of fixture placement angle up to 40°
 1.5mm lower vertical dimension
 Composition of multiple attachments with stable retention

Connect with exclusive driver
 Recommended tightening torque: 30Ncm



064 ZENEX
MULTI

Narrow DØ3.7

G/H (mm)	1.0	2.0	3.0	4.0	5.0	6.0
	G351	G352	G353	G354	G355	G356


M/R/W DØ3.7

G/H (mm)	1.0	2.0	3.0	4.0	5.0	6.0
	SQ351	SQ352	SQ353	SQ354	SQ355	SQ356

Kerator Angled Abutment

1.5mm lower vertical dimension
Composition of multiple attachments with stable retention

Connect with 1.2 Hex driver
Recommended tightening torque: 30Ncm



Narrow DØ3.7

G/H (mm) ——— 1.5 ——— 3.0



G35AN1.5



G35AN3.0

M/R/W DØ3.7

G/H (mm) ——— 1.5 ——— 3.0



SQ35AN1.5



SQ35AN3.0

Kerator CDPH Set

Composition: White spacer / denture cap connected black processing male & Replacement male (Red, Blue & Pink)

Select replacement male of adequate retention to use according to cases

Replace replacement male by using Kerator Magic tool



CDPH SET

066 ZENEX
09 MULTI

Kerator Male Cap

Retention power of KERATOR male cap is down (up to 20%) and red cap(angle) is included in the male package

Colors of cap determine retention power and it minimizes Denture Repair even there is any continuous bone loss



Kerator Magic Tool

Use to connect and remove replacement male on denture cap
As separated into two pieces, it is easy to maintenance



KMT002

067
ZENEX
MULTI

Kerator Torque Tip

Torque driver for Kerator Abutment

Round Type



KMD719

Kerator Impression Coping

Pick-Up Impression Coping for Kerator abutment
Closed Tray Type



DKI4845

Kerator Lab Analog

Lab Analog for Kerator Abutment



DKA3854

IZEN
I M P L A N T

069
ZENEX
MULTI

Prosthetic Flow Chart IV



Ti Link Abutment

72p




Ti Blank Abutment

73p



Fixture Lab Analog

47p



Pick-Up Impression Coping

43p



Transfer Impression Coping

45p



Cover Screw

25p



Healing Abutment

26p



1.2 Hex Driver

143p



ZENEX MULTI Fixture

22p


070 ZENEX
MULTI

ScanBody

ScanBody for manufacturing customized Titanium abutment
Use for Oral scan (Model scan available as well)

Tighten with 1.2 Hex Driver by hand

To order as a set as 'ScanBody + ScanBody Screw'
: Product Code + S (ex: ISBR4310S)



Narrow
DØ4.0

L : 14mm



TSBM4014



Mini
DØ4.0

L : 14mm



ISBM4014



R/W
DØ4.3

L : 10mm



ISBR4310




Ti Link Abutment

Abutment for manufacturing Cement/Combination-retained type prosthesis
 For manufacturing custom abutment (Titanium & Zirconia) and crown by CAD/CAM equipment
 Select specification fits for fixture Connection
 Use exclusive library for ZENEX MULTI Implant system
 Fixture Level Impression

Tighten with 1.2 Hex Driver
 Recommended tightening torque: 30Ncm

To order as a set as 'Abutment + Screw'
 : Product Code + S (ex: TLASRH4520S)



072 ZENEX MULTI

Narrow DØ4.0

G/H (mm)	1.0		2.0	
P/H (mm)	3.8			
	Hex	Non-Hex	Hex	Non-Hex
	TLATMH4010	TLATMN4010	TLATMH4020	TLATMN4020

M/R/W DØ4.5

G/H (mm)	1.0		2.0	
P/H (mm)	3.8			
	Hex	Non-Hex	Hex	Non-Hex
	TLASRH4510	TLASRN4510	TLASRH4520	TLASRN4520

Ti Blank Abutment


Manufacturing customized abutment with milling machine
 Select specification fits for fixture Connection
 Digital Impression

Tighten with 1.2 Hex Driver

Recommended tightening torque: 30Ncm

Product line-up applied for various milling machine brands
 (Milling machine manufacturer: Doowon, Manix, Vatech, RND)

To order as a set as 'Abutment + Screw'
 : Product Code + S (ex: TBASRH10AS)



073 ZENEX MULTI

For Arum Narrow DØ10



M/R/W DØ10



For Manix
Narrow
DØ10



M/R/W
DØ10



ZENEX KIT

COMPOSITION

134 ZENEX SURGERY KIT

150 ZENEX EASY SURGERY KIT

166 ZENEX SUPERWIDE KIT

176 ZENEX PROSTHETIC KIT


OTHER SURGICAL INSTRUMENTS

164 Tissue Punch

164 Trepine Drill

165 Bone Profiler


ZENEX SURGERY KIT



134 ZENEX
KIT

ZENEX SURGERY KIT

Composition on lower kit tray




135 ZENEX KIT

ZENEX SURGERY KIT Instruments

Guide Drill

- Drill to mark base hole on the bone for easy initial drilling
- Available to control drilling depth by using with stopper
- Product Code in the Kit: GD3F20L


D (mm)	Ø1.5 / Ø2.0	
Short	GD3F20S	
Long	GD3F20L	



Lindemann Drill

- Outward blades allow side cutting during drilling
- Available to use when removing ridge part of extraction socket
- Product Code in the Kit: LD4F25L

D (mm)	Ø1.5 / Ø2.0	Ø2.0 / Ø2.5
Short	LD4F20S	LD4F25S
Long	LD4F20L	LD4F25L



D (mm)	Ø3.0 / Ø3.5
Medium	LD4F35M

Drill Extension

- Lengthen surgical instrument in use by connecting to drill and Implant motor
- Be sure to use after accurate connection to prevent product detachment or bending that connected to Drill Extension


DE4513



Twist Drill

- Straight drill to shorten surgery steps
- Available to drill to desired depth by attaching stopper (marking line matches with stopper connection)
- Be sure to use stopper since depth control may be difficult due to excellent cutting power.
- Drill of 3.0F (Ø2.2/ 2.5) is applied only for X-Narrow of ZENEX PLUS I-System. (Product code in the Kit : TD2F25L)
- Drill of 3.0F (Ø2.2/ 2.7) is applied only for Narrow of ZENEX MULTI I-System. (Product code in the Kit : TD2F27L)

*Optional



137
ZENEX
KIT

D (mm)	2.2D (Ø1.8/ Ø2.2)	3.0F (Ø2.2 / Ø2.5)	3.0F (Ø2.2 / Ø2.7)	3.5F (Ø2.2 / Ø3.0)
Short	TD2F22S	TD2F25S*	TD2F27S*	TD2F30S
Long	TD2F22L	TD2F25L	TD2F27L	TD2F30L
D (mm)	4.0F (Ø3.0 / Ø3.6)	4.5F (Ø3.0 / Ø3.6 / Ø4.1)	5.0F (Ø3.0 / Ø4.1 / Ø4.6)	
Short	TD2F36S	TD2F41S	TD2F46S	
Long	TD2F36L	TD2F41L	TD2F46L	

ZENEX SURGERY KIT Instruments

Step Drill


- Straight drill to shorten surgery steps
- As stopper is applied in the drill, secure to use although depth control may be difficult due to excellent cutting power
- Colored line on grip part of drill indicates diameter of drill and matching fixture size

- Two kinds of 5.5F Drill (Ø 4.6 / Ø 5.2) are included.

5.5F
L (mm) (Ø4.6 / Ø5.2)

11.5 TD3F5211

15.0 TD3F5215



- Optional extras: Three kinds of step drill for 20mm extra long



3.5F
L (mm) (Ø2.2 / Ø3.0)

20.0 TD2F3020



4.0F
L (mm) (Ø3.0 / Ø3.6)

TD2F3620




4.5F
L (mm) (Ø3.0 / Ø3.6 / Ø4.1)

TD2F4120

Drill Stopper


- Indicated number on each stopper means exposed length, excluding head of drill (Y-Dim) when connected to drill
- Differentiate color by length for easy identification



139 ZENEX KIT

Cortical Drill





- Drill used for removing cortical bone from hard bone
- Drill to the bottom of each marking line
- The lower marking line is for drilling on normal bone and the upper marking line is for drilling on hard bone



ZENEX SURGERY KIT Instruments

Direction Pin

- Confirm gingiva height and path after implant placement
- DPSR4509: Available to use in case of ZENEX I-System (Regular & Wide) and ZENEX T-System (Regular & Wide)
- *Only use red marked components for I-System**


I-System	X-Narrow (Ø3.5)	Narrow (Ø3.5)	Mini (Ø3.5)	R/W (Ø4.5)
				
	DPSX3510*	DPSN3510*	DPSM3510*	DPSR4509*
COLOR	Red	Purple	Green	Metal

- For ZENEX MULTI I-System (based on connection): Narrow, Mini, Regular & Wide
- For ZENEX PLUS I-System (based on connection): X-Narrow, Mini, Regular & Wide

T-System	Mini (Ø3.5)	R/W (Ø4.5)	R-System	M/R/W (Ø4.0)
				
	DPTM4010	DPSR4509		DPRR4009
COLOR	Yellow	Metal	COLOR	Blue

Parallel Pin

- Used for confirmation of drilling depth and path
- Available to use after drilling with Ø2.2 Drill and Ø3.0 drill
- Composition in the kit: PP352230 (1 ea) + PP402230 (2 ea)



D (mm)	Ø3.5	Ø4.0	Ø5.0	Ø6.0
	PP352230	PP402230	PP502230	PP602230

Fixture Driver

- For Contra Angle: Directly attached to fixture by connecting to handpiece when using Implant motor
- For Wrench: Directly attached to fixture by using Torque Wrench
- Two-staged coloring enables measurement up to 5mm height at 1mm intervals
- Be sure to use after accurate connection because once used with incorrectly connected driver, it is impossible to remove fixture.
- Be sure to use accurate fixture driver according to the system of fixture.
- ※ Fixture Driver for Narrow of ZENEX MULTI I-System: Same with Mini of ZENEX T-System (Optional)
- ※ Code not included in the FDCSX17S, FDRSX17S, All four drivers for Narrow of ZENEX MULTI I-System (FDCTM21S, FDCTM21L, FDRTM21S, FDRTM21L)
- *Only use red marked components for I-System

I-System

Contra Angle



X-Narrow Narrow M/R/W

Wrench



X-Narrow Narrow M/R/W

	X-Narrow	Narrow	M/R/W	X-Narrow	Narrow	M/R/W
Short	FDCSX17S*	FDCTM21S	FDCSR25S*	FDRSX17S*	FDRTM21S	FDRSR25S*
Long	FDCSX17L*	FDCTM21L*	FDCSR25L*	FDRSX17L*	FDRTM21L*	FDRSR25L*
COLOR	Red	Yellow	Green	Red	Yellow	Green

- For ZENEX MULTI I-System (based on connection): Narrow, Mini, Regular & Wide
- For ZENEX PLUS I-System (based on connection): X-Narrow, Mini, Regular & Wide

ZENEX SURGERY KIT Instruments

T-System

Contra Angle



Mini

R/W

Wrench



Mini

R/W

Short

Long

COLOR

FDCTM21S

FDCTR25S

FDRTM21S

FDRTR25S

FDCTM21L

FDCTR25L

FDRTM21L

FDRTR25L

Yellow

Black

Yellow

Black

142 ZENEX KIT

R-System

Contra Angle



M/R/W

Wrench



M/R/W

Short

Long

COLOR

FDCRR23S

FDRRR23S

FDCRR23L


FDRRR23L

Blue

Blue

1.2 Hex Driver

- Driver for connecting with Torque Wrench
- Used for connecting to cover screw and surgical instrument with 1.2 Hex
- Follow recommended tightening torque calculation value
- Be sure to check transformation of Hex tip before use



Torque Wrench

- Used for placing fixture and tightening implant abutments and screws
- Confirm torque calculation value by pulling the bar to be aligned with the marked line
- Recommended to use by pressing head to prevent detachment of Torque Wrench and surgical instruments



TW60B

ZENEX SURGERY KIT Instruments

Depth Gauge

- Measuring drilling depth
- Measurable depth: 6~20mm



DG1720

Steel Bowl

- Used for cleaning surgery instruments




STBOWL

Drilling Sequence

Ø 3.0mm Fixture

(ZENEX PLUS I-System X-Narrow)




Bone Quality	Guide Drill	Twist Drill (Ø2.2)	Twist Drill (Ø2.5)	Cortical Drill	Ø3.0 Fixture
Soft	▶	▶			Implant Placement
Normal	▶	▶	▶		
Hard	▶	▶	▶	▶	

145
ZENEX
KIT

Ø 3.0mm Fixture

(ZENEX MULTI I-System Narrow)




Bone Quality	Guide Drill	Twist Drill (Ø2.2)	Twist Drill (Ø2.7)	Cortical Drill	Ø3.0 Fixture
Soft	▶	▶			Implant Placement
Normal	▶	▶	▶		
Hard	▶	▶	▶	▶	

Drilling Sequence

Ø 3.5mm Fixture


(ZENEX MULTI & PLUS I-System Mini)



Bone Quality	Guide Drill	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Cortical Drill (Lower Line)	Cortical Drill (Upper Line)	Ø3.5 Fixture
Soft	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶		
Hard	▶	▶	▶		▶	

Ø 4.0mm Fixture

(ZENEX MULTI & PLUS I-System Regular)




Bone Quality	Guide Drill	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Twist Drill (Ø3.6)	Cortical Drill (Lower Line)	Cortical Drill (Upper Line)	Ø4.0 Fixture
Soft	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶		▶	

Drilling Sequence

Ø 4.5mm Fixture

(ZENEX MULTI & PLUS I-System Regular)




Bone Quality	Guide Drill	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Twist Drill (Ø3.6)	Twist Drill (Ø4.1)	Cortical Drill (Lower Line)	Cortical Drill (Upper Line)	Ø4.5 Fixture
Soft	▶	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶		▶	

147
ZENEX
KIT

Ø 5.0mm Fixture

(ZENEX MULTI & PLUS I-System Wide)




Bone Quality	Guide Drill	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Twist Drill (Ø4.1)	Twist Drill (Ø4.6)	Cortical Drill (Lower Line)	Cortical Drill (Upper Line)	Ø5.0 Fixture
Soft	▶	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶		▶	

Drilling Sequence

Ø 5.5mm Fixture

(ZENEX MULTI & PLUS I-System Wide)




Bone Quality	Guide Drill	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Twist Drill (Ø4.6)	Step Drill (Ø5.2)	Cortical Drill (Lower Line)	Cortical Drill (Upper Line)	Ø5.5 Fixture
Soft	▶	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶		▶	

IZEN
I M P L A N T

149
ZENEX
KIT


ZENEX EASY SURGERY KIT



150 ZENEX
KIT

ZENEX EASY SURGERY KIT

Composition on lower kit tray



151 ZENEX KIT

ZENEX EASY SURGERY KIT Instruments

Guide Drill

- Drill to mark base hole on the bone for easy initial drilling
- Available to control drilling depth by using with stopper
- Product Code in the Kit: GD3F20L

D (mm) Ø1.5 / Ø2.0

Short

GD3F20S

Long

GD3F20L



Lindemann Drill

- Outward blades allow side cutting during drilling
- Available to use when removing ridge part of extraction socket
- Product Code in the Kit: LD4F25L

D (mm) Ø1.5 / Ø2.0 Ø2.0 / Ø2.5

Short

LD4F20S

LD4F25S

Long


LD4F20L

LD4F25L

D (mm) Ø3.0 / Ø3.5

Medium

LD4F35M



Drill Extension

- Lengthen surgical instrument in use by connecting to drill and Implant motor
- Be sure to use after accurate connection to prevent product detachment or bending that connected to Drill Extension


DE4513



Step Drill

- Straight drill to shorten surgery steps
- As stopper is applied in the drill, secure to use although depth control may be difficult due to excellent cutting power
- Colored line on grip part of drill indicates diameter of drill and matching fixture size.

* Optional



153 ZENEX KIT

- Two kinds of step drill for X-Narrow (ZENEX PLUS I-System) & Narrow (ZENEX MULTI I-System) are included in the Kit.



X-Narrow

Narrow

3.0F (Ø2.2 / Ø2.5)

3.0F (Ø2.2 / Ø2.7)

L (mm)

15.0

TD2F2515

TD2F2715

ZENEX EASY SURGERY KIT Instruments

Step Drill

* Optional

L (mm)	2.2D	3.5F	4.0F
	(Ø1.8 / Ø2.2)	(Ø2.2 / Ø3.0)	(Ø3.0 / Ø3.6)
6.0	TD2F2206 *	TD2F3006 *	TD2F3606 *
7.0	TD2F2207	TD2F3007	TD2F3607
8.5	TD2F2208	TD2F3008	TD2F3608
10.0	TD2F2210	TD2F3010	TD2F3610
11.5	TD2F2211	TD2F3011	TD2F3611
13.0	TD2F2213	TD2F3013	TD2F3613
15.0	TD2F2215 *	TD2F3015 *	TD2F3615 *
18.0	-	TD2F3018 *	TD2F3618 *
20.0	-	TD2F3020 *	TD2F3620 *
L (mm)	4.5F	5.0F	5.5F
	(Ø3.0 / Ø3.6 / Ø4.1)	(Ø3.0 / Ø4.1 / Ø4.6)	(Ø4.6 / Ø5.2)
6.0	TD2F4106 *	TD2F4606 *	TD3F5206 *
7.0	TD2F4107	TD2F4607	TD3F5207
8.5	TD2F4108	TD2F4608	TD3F5208
10.0	TD2F4110	TD2F4610	TD3F5210
11.5	TD2F4111	TD2F4611	TD3F5211
13.0	TD2F4113	TD2F4613	TD3F5213
15.0	TD2F4115 *	TD2F4615 *	TD3F5215 *
18.0	TD2F4118 *	-	-
20.0	TD2F4120 *	-	-

Cortical Drill

- Drill used for removing cortical bone from hard bone
- Drill to the bottom of each marking line
- The lower marking line is for drilling on normal bone and the upper marking line is for drilling on hard bone



	3.0C (X-Narrow)	3.0C (Narrow)	3.5C	4.0C	4.5C	5.0C	5.5C
	CD3FXF30	CD3FNF30	CD4FF35	CD4FF40	CD4FF45	CD4FF50	CD4FF55
COLOR	None	Purple	Yellow	Green	Blue	Red	Yellow

Direction Pin

- Confirm gingiva height and path after implant placement
- DPSR4509: Available to use in case of ZENEX I-System (Regular & Wide) and ZENEX T-System (Regular & Wide)
- *Only use red marked components for I-System**

I-System	X-Narrow (Ø3.5)	Narrow (Ø3.5)	Mini (Ø3.5)	R/W (Ø4.5)
COLOR	DPSX3510* Red	DPSN3510* Purple	DPSM3510* Green	DPSR4509* Metal


- For ZENEX MULTI I-System (based on connection): Narrow, Mini, Regular & Wide
- For ZENEX PLUS I-System (based on connection): X-Narrow, Mini, Regular & Wide

T-System	Mini (Ø3.5)	R/W (Ø4.5)	R-System	M/R/W (Ø4.0)
COLOR	DPTM4010 Yellow	DPSR4509 Metal	COLOR	DPRR4009 Blue

ZENEX EASY SURGERY KIT Instruments

Parallel Pin

- Used for confirmation of drilling depth and path
- Available to use after drilling with Ø2.2 Drill and Ø3.0 drill
- Composition in the kit: PP352230 (1 ea) + PP402230 (2 ea)



D (mm)	Ø3.5	Ø4.0	Ø5.0	Ø6.0
	PP352230	PP402230	PP502230	PP602230

Fixture Driver

- For Contra Angle: Directly attached to fixture by connecting to handpiece when using Implant motor
- For Wrench: Directly attached to fixture by using Torque Wrench
- Two-staged coloring enables measurement up to 5mm height at 1mm interval
- Be sure to use after accurate connection because once used with incorrectly connected driver, it is impossible to remove fixture.
- Be sure to use accurate fixture driver according to the system of fixture.
- ※ Fixture Driver for Narrow of ZENEX MULTI I-System: Same with Mini of ZENEX T-System (Optional)
- ※ Code not included in the FDCSX17S, FDRSX17S, All four drivers for Narrow of ZENEX MULTI I-System (FDCTM21S, FDCTM21L, FDRTM21S, FDRTM21L)
- *Only use red marked components for I-System

I-System

Contra Angle



X-Narrow Narrow M/R/W

Wrench



X-Narrow Narrow M/R/W

Short

FDCSX17S* FDCTM21S FDCSR25S*

FDRSX17S* FDRTM21S FDRSR25S*

Long

FDCSX17L* FDCTM21L* FDCSR25L*

FDRSX17L* FDRTM21L* FDRSR25L*

COLOR

Red

Yellow

Green

Red

Yellow

Green

- For ZENEX MULTI I-System (based on connection): Narrow, Mini, Regular & Wide
- For ZENEX PLUS I-System (based on connection): X-Narrow, Mini, Regular & Wide

T-System

Contra Angle



Mini

R/W

Wrench



Mini

R/W

Short

Long

COLOR

FDCTM21S

FDCTR25S

FDRTM21S

FDRTR25S

FDCTM21L

FDCTR25L

FDRTM21L

FDRTR25S

Yellow

Black

Yellow

Black

R-System

Contra Angle



M/R/W

Wrench



M/R/W

Short

Long

COLOR

FDCRR23S

FDRRR23S

FDCRR23L

FDRRR23L


Blue

Blue

1.2 Hex Driver

- Driver for connecting with Torque Wrench
- Used for connecting to cover screw and surgical instrument with 1.2 Hex
- Follow recommended tightening torque calculation value
- Be sure to check transformation of Hex tip before use

L (mm)	10.0	15.0
	HEDR1210	HEDR1215



Torque Wrench

- Used for placing fixture and tightening implant abutments and screws
- Confirm torque calculation value by pulling the bar to be aligned with the marked line
- Recommended to use by pressing head to prevent detachment of Torque Wrench and surgical instruments



TW60B

Depth Gauge

- Measuring drilling depth
- Measurable depth: 6~20mm



DG1720

Steel Bowl

- Used for cleaning surgery instruments




STBOWL

Drilling Sequence

Ø 3.0mm Fixture


(ZENEX PLUS I-System X-Narrow)



Bone Quality	Guide Drill	Step Drill (Ø2.2)	Step Drill (Ø2.5)	Cortical Drill	Ø3.0 Fixture
Soft	▶	▶	▶		Implant Placement
Normal	▶	▶	▶		
Hard	▶	▶	▶	▶	

Ø 3.0mm Fixture

(ZENEX MULTI I-System Narrow)




Bone Quality	Guide Drill	Step Drill (Ø2.2)	Step Drill (Ø2.7)	Cortical Drill	Ø3.0 Fixture
Soft	▶	▶	▶		Implant Placement
Normal	▶	▶	▶		
Hard	▶	▶	▶	▶	

Drilling Sequence

Ø 3.5mm Fixture

(ZENEX MULTI & PLUS I-System Mini)




Bone Quality	Guide Drill	Step Drill (Ø2.2)	Step Drill (Ø3.0)	Cortical Drill (Lower Line)	Cortical Drill (Upper Line)	Ø3.5 Fixture
Soft	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶		
Hard	▶	▶	▶		▶	

161 ZENEX KIT

Ø 4.0mm Fixture

(ZENEX MULTI & PLUS I-System Regular)




Bone Quality	Guide Drill	Step Drill (Ø2.2)	Step Drill (Ø3.0)	Step Drill (Ø3.6)	Cortical Drill (Lower Line)	Cortical Drill (Upper Line)	Ø4.0 Fixture
Soft	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶		▶	

Drilling Sequence

Ø 4.5mm Fixture


(ZENEX MULTI & PLUS I-System Regular)



Bone Quality	Guide Drill	Step Drill (Ø2.2)	Step Drill (Ø3.0)	Step Drill (Ø3.6)	Step Drill (Ø4.1)	Cortical Drill (Lower Line)	Cortical Drill (Upper Line)	Ø4.5 Fixture
Soft	▶	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶		▶	

Ø 5.0mm Fixture

(ZENEX MULTI & PLUS I-system Wide)




Bone Quality	Guide Drill	Step Drill (Ø2.2)	Step Drill (Ø3.0)	Step Drill (Ø4.1)	Step Drill (Ø4.6)	Cortical Drill (Lower Line)	Cortical Drill (Upper Line)	Ø5.0 Fixture
Soft	▶	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶		▶	

Drilling Sequence

Ø 5.5mm Fixture

(ZENEX MULTI & PLUS I-System Wide)




Bone Quality	Guide Drill	Step Drill (Ø2.2)	Step Drill (Ø3.0)	Step Drill (Ø4.6)	Step Drill (Ø5.2)	Cortical Drill (Lower Line)	Cortical Drill (Upper Line)	Ø5.5 Fixture
Soft	▶	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶		▶	

Other Surgical Instruments

Tissue Punch

- Marked with a gap of 2mm each to measure gingiva height
- Surgical instrument for Flapless surgery



D (mm)	Ø3.3	Ø3.8	Ø4.3	Ø4.8	Ø5.3
	TP3308	TP3808	TP4308	TP4808	TP5308

164 ZENEX KIT

Trephine Drill

- Use for collecting autogenous bone or removing failed fixture
- In case of Superwide fixture placement, available to use as initial drill



D (mm)	Ø3.7/4.5	Ø4.2/5.0	Ø4.7/5.5	Ø5.2/6.0	Ø5.7/6.5	Ø6.2/7.0
	TD3718	TD4218	TD4718	TD5218	TD5718	TD6218

Bone Profiler

- Use for removing marginal bone around fixture during the one-stage and two-stages surgery- Available to connect Healing Abutment by removing marginal bone after connecting Bone Profiler Guide Screw to the connection on Fixture
- Secure guided length of Bone Profiler Guide Screw to remove bone accurately for deeply placed fixture
- Bone Profiler Guide Screw is used for protecting fixture connection during use



I-System ※ Only for I-System

D (mm)	Ø4.1 (For Ø3.0 X-Narrow)		Ø4.1 (For Ø3.0 Narrow)	
	SBP41SXS		SBP41SNS	
D (mm)	Ø4.1	Ø4.6	Ø5.6	Ø7.6
	SBP41SS	SBP46SS	SBP56SS	SBP76SS

T-System ※ Only for T-System

D (mm)	Ø4.1	Ø4.6	Ø5.6	Ø7.6
	SBP41TS	SBP46TS	SBP56TS	SBP76TS

R-System ※ Only for R-System

D (mm)	Ø4.1	Ø4.6	Ø5.6	Ø7.6
	SBP41RS	SBP46RS	SBP56RS	SBP76RS

	I-System		T-System		R-System		
Bone Profiler Guide Screw	X-Narrow	Narrow	Mini	Regular	Mini	Mini	Regular
	SBPGSSX14	SBPGSSN20	SBPGSSM20	SBPGSSR20	SBPGSTM20	SBPGSRM18	SBPGSRR18


ZENEX SUPERWIDE KIT

166 ZENEX
KIT



ZENEX SUPERWIDE KIT

Composition on lower kit tray



167
ZENEX
KIT

ZENEX SUPERWIDE KIT Instruments


Guide Drill

- Drill to mark base hole on the bone for easy initial drilling
- Available to control drilling depth by using with stopper
- Product Code in the Kit: GD3F20L

D (mm) $\overline{\quad\quad\quad}$ $\text{\O}1.5 / \text{\O}2.0$

Short GD3F20S

Long GD3F20L



168 ZENEX KIT

Lindemann Drill

- Outward blades allow side cutting during drilling
- Available to use when removing ridge part of extraction socket
- Product Code in the Kit: LD4F25L


D (mm) $\overline{\quad\quad\quad}$ $\text{\O}1.5 / \text{\O}2.0$ $\text{\O}2.0 / \text{\O}2.5$

Short LD4F20S LD4F25S

Long LD4F20L LD4F25L

D (mm) $\overline{\quad\quad\quad}$ $\text{\O}3.0 / \text{\O}3.5$

Medium LD4F35M



Drill Extension


- Lengthen surgical instrument in use by connecting to drill and Implant motor
- Be sure to use after accurate connection to prevent product detachment or bending that connected to Drill Extension

DE4513



Step Drill

- Straight drill to shorten surgery steps
- As stopper is applied in the drill, secure to use although depth control may be difficult due to excellent cutting power
- Colored line on grip part of drill indicates diameter of drill and matching fixture size.
- The length of 13mm step drill is applied only in the kit.
- ※ TD3F6213: Only use for drilling on soft bone when Ø 7.0 fixture is placed.



169
ZENEX
KIT

D (mm)	2.2D (Ø1.8 / Ø2.2)	3.5F (Ø2.2 / Ø3.0)	5.0F (Ø3.0 / Ø4.1 / Ø4.6)
	TD2F2213	TD2F3013	TD2F4613
D (mm)	5.5F (Ø4.6 / Ø5.2)	6.0F (Ø4.6 / Ø5.5)	7.0F (Ø5.5 / Ø6.2)
	TD3F5213	TD3F5513	TD3F6213
D (mm)	7.0F (Ø5.5 / Ø6.5)		
	TD3F6513		

ZENEX SUPERWIDE KIT Instruments

Cortical Drill

- Drill used for removing cortical bone from hard bone
- Drill to the bottom of each marking line
- In case of Wide fixture (more than Ø6.0) placement, use only for hard bone case.




Direction Pin

- Confirm gingiva height and path after implant placement
 - DPSR4509: Available to use in case of ZENEX I-System (Regular & Wide) and ZENEX T-System (Regular & Wide)
 - 2 pins of each system are included in the kit, in accordance with fixture system.
- The user may choose the product, matching with the system



Parallel Pin

- Used for confirmation of drilling depth and path
- Available to use after drilling with Ø2.2 Drill and Ø3.0 drill
- Composition in the kit: PP402230 (2 ea)



D (mm)	Ø3.5	Ø4.0	Ø5.0	Ø6.0
	PP352230	PP402230	PP502230	PP602230

Fixture Driver

- For Contra Angle: Directly attached to fixture by connecting to handpiece when using Implant motor
- For Wrench: Directly attached to fixture by using Torque Wrench
- Two-staged coloring enables measurement up to 5mm height at 1mm interval
- Be sure to use after accurate connection because once used with incorrectly connected driver, it is impossible to remove fixture
- Be sure to use accurate fixture driver according to the system of fixture.
- Short type of fixture driver is included only

I-System * Only for I-System

Contra Angle

Wrench



M/R/W

M/R/W

Short
COLOR

FDCSR25S
Green

FDRSR25S
Green

ZENEX SUPERWIDE KIT Instruments

T-System ※ Only for T-System

Contra Angle



R/W

Short
COLOR

FDCTR25S
Black

Wrench



R/W

FDRTR25S
Black

R-System ※ Only for R-System

Contra Angle



M/R/W

Short
COLOR

FDCRR23S
Blue

Wrench




M/R/W

FDRRR23S
Blue

1.2 Hex Driver

- Driver for connecting with Torque Wrench
- Used for connecting to cover screw and surgical instrument with 1.2 Hex
- Follow recommended tightening torque calculation value
- Be sure to check transformation of Hex tip before use



Torque Wrench

- Used for placing fixture and tightening implant abutments and screws
- Confirm torque calculation value by pulling the bar to be aligned with the marked line
- Recommended to use by pressing head to prevent detachment of Torque Wrench and surgical instruments



TW60B

ZENEX SUPERWIDE KIT Instruments

Depth Gauge

- Measuring drilling depth
- Measurable depth: 6~20mm




DG1720

Drilling Sequence

Ø 6.0mm Fixture

(ZENEX MULTI I-System Wide)




Bone Quality	Guide Drill	Step Drill (Ø2.2)	Step Drill (Ø3.0)	Step Drill (Ø4.6)	Step Drill (Ø5.2)	Step Drill (Ø5.5)	Cortical Drill	Ø6.0 Fixture
Soft	▶	▶	▶	▶	▶ (6mm)			Implant Placement
Normal	▶	▶	▶	▶		▶ (6mm)		
Hard	▶	▶	▶	▶		▶ (6mm)	▶	

175 ZENEX KIT

Ø 7.0mm Fixture

(ZENEX MULTI I-System Wide)



Bone Quality	Guide Drill	Step Drill (Ø2.2)	Step Drill (Ø3.0)	Step Drill (Ø4.6)	Step Drill (Ø5.5)	Step Drill (Ø6.2)	Step Drill (Ø6.5)	Cortical Drill	Ø7.0 Fixture
Soft	▶	▶	▶	▶	▶	▶ (6mm)			Implant Placement
Normal	▶	▶	▶	▶	▶		▶ (6mm)		
Hard	▶	▶	▶	▶	▶		▶ (6mm)	▶	

ZENEX PROSTHETIC KIT

176 ZENEX
KIT




ZENEX PROSTHETIC KIT

Composition on lower kit tray



Kerator Magic Tool
KMT002




177
ZENEX
KIT

1.2 Hex Driver

- Driver for connecting with Torque Wrench
- Used for connecting to cover screw and surgical instrument with 1.2 Hex
- Follow recommended tightening torque calculation value
- Be sure to check transformation of Hex tip before use


L (mm)	10.0	15.0	25.0
	HEDR1210	HEDR1215	HEDR1225



1.2 Hex Machine Screw Driver

- Driver for implant motor
- Implement holding instrument
- Available to use by hand when connected to Finger Driver for Driver Shaft

L (mm)	10.0	15.0
	HEDM1210	HEDM1215



Abutment Remover

- Use for removing abutment when it is impossible to remove abutment by hands due to strong connection force.
 - Remove Abutment by keep tightening the remover clockwise after removing abutment screw
- ※ R-System use only



ARRM18

Finger Driver for Driver Shaft

- Available to rotate with hand by connecting to every surgery instrument for motor



FDDS100

ZENEX PROSTHETIC KIT Instruments

Ball Abutment Driver

- Torque driver for Ball Abutment



BAD24

Multi Straight Abutment Driver

- Torque driver for Multi Straight Abutment



MSADSR20

Torque Wrench

- Used for placing fixture and tightening implant abutments and screws
- Confirm torque calculation value by pulling the bar to be aligned with the marked line
- Recommended to use by pressing head to prevent detachment of Torque Wrench and surgical instruments



TW60B

Kerator Magic Tool

- Use to connect and remove replacement male on denture cap
- As separated into two pieces, it is easy for maintenance (Optional)



KMT002

Kerator Torque Tip

- Torque driver for Kerator Abutment (Optional)

Round Type



KMD719

IZEN
I M P L A N T

183
ZENEX
KIT

1. Description of Izenimplant System's product

This product is a dental implant which is put into the alveolar bone in order to support, or maintain the prosthetic tooth or denture when a patient's teeth are partially or totally lost. To enhance the osseointegration with the alveolar bone, this titanium dental implant is treated with SLA surface.

2. Indications for use

The ZENEX MULTI implant and ZENEX PLUS implant are indicated for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple unit restorations. It is intended for delayed loading. Products with diameter equal to or above 5mm and length below 7mm is intended to be used in the molar region.

3. How to store

Keep it in a dry place having no direct sunlight and room temperature (1~30°C)

4. How to use

4.1 Preliminary preparations

- ① For a successful dental implant treatment, careful evaluation should be made for each patient from the medical or dental perspective and clinical and radiologic analysis and model analysis are required.
- ② The internal package of the dental implant should be checked. The internal package of the dental implant is sterilized before delivered. So, if any seal is torn or opened, the product should not be used.

4.2 How to operate

【1st surgery】

- ① Make an incision in the soft tissue and drill it in a way that a fixture can be put into it. (Drilling : 1,200~1,500rpm, Cortical Drilling : 500~800rpm)

Ø3.0 F	GD→Ø2.2D→Ø2.5D or 2.7D
Ø3.5 F	GD→Ø2.2D→Ø3.0D→CD
Ø4.0 F	GD→Ø2.2D→Ø3.0D→Ø3.6D→CD
Ø4.5 F	GD→Ø2.2D→Ø3.0D→Ø3.6D→Ø4.1D→CD
Ø5.0 F	GD→Ø2.2D→Ø3.0D→Ø4.1D→Ø4.6D→CD
Ø5.5 F	GD→Ø2.2D→Ø3.0D→Ø4.6D→Ø5.2SD→CD
Ø6.0 F	GD→Ø2.2SD→Ø3.0SD→Ø4.6SD→Ø5.2SD→Ø5.5SD
Ø7.0 F	GD→Ø2.2SD→Ø3.0SD→Ø4.6SD→Ø5.5SD→Ø6.2SD→Ø6.5SD

GD: Guide Drill, D:Twist/Step Drill, SD:Step Drill, CD:Cortical Drill

*Remark : Drilling is omissible depending on extracted socket size, if bone density is hard, it is recommended to use Cortical Drill properly.

- ② Check out the type and size of Fixture specified in the package label and take out the double packaged ampoule.
- ③ Turn the ampoule cap by hand to open and find the Fixture which is fixed to the titanium bushing. Connect the Fixture Driver to the Handpiece or Ratchet.
- ④ Take the Fixture to the mouth and put it into the drilled hole while making sure that any metal or saliva does not contaminate it. While changing the torque value of the Handpiece from low value to high value, put the Fixture into the alveolar bone in 20rpm for full implantation (The recommended torque value for implantation : 30 ~ 45Ncm)
- ⑤ Use the Hex Driver to fix the Cover Screw or Healing Abutment to the Fixture with 5~8Ncm before making a suture in the surgery area.

【2nd surgery】

- ① Once the alveolar mucosa is healed and the osseointegration is made, make an incision in the upper soft tissue above the implanted area to expose the Cover Screw.
- ② Remove the Cover Screw and connect the Healing Abutment to the Fixture.
- ③ Suture the soft tissue around the Healing Abutment.

4.3 How to store and control after use

- ① As all products put into the mouth are intended for disposable use, they are not allowed for reuse.
- ② Once opened, any product should be disposed of.
- ③ Any packing materials for the product used should be disposed of.
- ④ The surgical tool which is used should be immediately cleaned by using alcohol, cleaning agent or distilled water and then sterilized and dried before stored.

5. Warning

- ① The selection of an improper patient and a wrong surgery may lead to a failure of the treatment.
- ② As the product is a disposable one, it should not be reused.
- ③ Risks from reusing the disposable medical device
 - ▶ The reuse has not been verified so it may cause an infection.
 - ▶ The product's defective connection, jammed abutment, deformation or screw loosening may occur.
- ④ The manufacturer does not take any responsibility for the product which is re-sterilized by the user.
- ⑤ The Izen implant System has not been evaluated for safety and compatibility in the MR (magnetic resonance) environment. The Izen implant System has not been tested for heating or migration, or image artifact in the MR environment. The safety of fixture in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

6. Prohibitions

The treatment should not be done for the following patients.

- ① Patients having problems in coagulation or disability in the healing of bone or wound
- ② Patients having any uncontrolled diabetes, heavy smokers or alcoholics
- ③ Patients with decreased immune system due to the chemotherapy and radiotherapy.
- ④ Patients having oral infection or inflammation (or patients having improper oral hygiene or bruxism)
- ⑤ Patients having untreatable disorder in the occlusion or joint or having insufficient arch space
- ⑥ Other patients who are not suitable for treatment
- ⑦ Patients who are allergic or too sensitive to the raw material used for implant
- ⑧ Patients who are not cooperative in treatment

7. Adverse reactions (side effects)

- ① Bone loss, low primary stability, dental implant loss, inflammation or neural damage may occur, thus leading to a failure in the dental implant surgery.
- ② Local complications such as swelling, pain, phonetic difficulties, hematoma, hemorrhage, infection, inflammation, ulcer or opening of sutured area may occur.
- ③ Edema, separation of sutured area
- ④ Numbness from nerve damage after surgery is usually a temporary symptom, but some cases are permanent.
- ⑤ Permanent paresthesia, dysesthesia, loss of maxillary/mandibular ridge bone.
- ⑥ Oroantral or oronasal fistulae, unfavorably affected adjacent teeth, irreversible damage to adjacent teeth, fractures of implant, jaw, bone or prosthesis, aesthetic problems, exfoliation, hyperplasia.

8. General precautions

- ① Only a dentist who has received implant surgery training and practice should be allowed to use the product.
- ② The operator should fully understand how to use the product and other cautions before selecting the product which is proper for the treatment plan.
- ③ Patients who have excessive masticatory force may experience failed osseointegration, or loosening or breaking off of the dental implant.
So, make sure to perform the surgery with the sufficient number, thickness and length of dental implants.
- ④ Prior to the surgery, check out the surgery tool for abrasion.
- ⑤ Make sure to supply sufficient water during the surgery so that the bone does not suffer necrosis by the heat.
- ⑥ As the product is sterilized, make sure to check the product for any packing status, expiry date or damage.
- ⑦ As the product is sterilized, make sure to open it in a clean environment immediately prior to the surgery.
- ⑧ The selection and treatment of an improper patient may lead to a failure of the dental implant or bone loss around the implant area after the surgery.
- ⑨ If a patient has any bone disease (such as osteoporosis or osteomalacia) or bone metabolism disorder, consider this carefully prior to the treatment.
- ⑩ As an implant having a diameter of 4.25mm or below may get broken due to limited strength, its use is not recommended for the oral posterior area.
- ⑪ The use of an implant having a diameter of not less than 5mm and length below 7mm is limited to the posterior area and the dental implant shall be splinted with other prostheses.
- ⑫ The product having the diameter of 3.25mm or below shall be used only for the mandibular anterior region in order to avoid its fracture due to excessive load force.

“Disposable, no reuse, medical device”

※ For detailed information of each product, please refer to the catalog or our website.



 **Manufacturer : Izenimplant Co., Ltd.**

1, 2Dong, 26-32, Suworam 4-gil, Seotan-myeon, Pyeongtaek-si, Gyeonggi-do, 17703, Republic of Korea
Tel.82-31-662-0657 Fax.82-31-662-0658

Rx only
PRESCRIPTION ONLY

For USA Only : Federal law restricts this device to sale by or on the order of a dentist