

ZENEX  
IMPLANT SYSTEM

# GUIDE KIT

The Best way for the successful surgery



Fixture  
Size




ZENEX  
Sleeve  
Size



# ZENEX GUIDE KIT

ZENEX Implant Guide Kit for  $\Phi 3.5/\Phi 4.0/\Phi 4.5/\Phi 5.0$   
For Use with **Regular Size Sleeve**


**Off-Set** 9 / 10.5 / 12mm  
**Code** ZENGIKIT



# ZENEX GUIDE MINI KIT

ZENEX Implant Guide Kit for  $\Phi 3.0/\Phi 3.5$   
For Use with **Mini Size Sleeve**

**Off-Set** 12mm  
**Code** ZENGMKIT





# Fixture Direction

Control the fixture direction by matching the groove of the sleeve and the shape of the driver

